

Exercise 1: Spatial filters

Artifacts are a major issue in BCI. One way to reduce the impact of artifacts are spatial filters. In this exercise you will study the effects of different spatial filters on artifact-contaminated EEG.

The file `artefacts.mat` contains a segment of a 120 channel EEG recording that is contaminated with artifacts. Use the MATLAB function `load` to load the `.mat` file. The EEG is sampled at a rate of 500 Hz and band pass filtered between 1 Hz and 200 Hz. Electrode setup and channel information are shown on the next pages. Please compute the following spatial filters for positions Fz, Cz and Oz. Note: in the case of laplacian derivation, anterior means to the front, posterior means to the back, dexter means to the right and sinister to the left with respect to a chosen electrode position on which the spatial filter is applied.

Exercise 1.1: Bipolar derivation $X_{bip} = X_{anterior} - X_{posterior}$

Exercise 1.2: Laplace derivation $X_{lap} = X - \frac{1}{4} \left(X_{anterior} + X_{dexter} + X_{posterior} + X_{sinister} \right)$

Exercise 1.3: Common average reference derivation (CAR) $X_{CAR} = X - average(all\ channels)$


Task:

- Describe the effects of the different spatial filters on the EEG signal, especially on artifacts.

Exercise 2: Artifacts in EEG


Below are some examples of EEG independent components. Label them correctly and briefly explain your reasons. Note that some components might not be artefactual, but brain related.

1.




☐ Brain ☐ Muscle ☐ Eye ☐ Heart ☐ Line Noise ☐ Chan Noise ☐ Other ☐ ?

2.




☐ Brain ☐ Muscle ☐ Eye ☐ Heart ☐ Line Noise ☐ Chan Noise ☐ Other ☐ ?

3.




☐ Brain ☐ Muscle ☐ Eye ☐ Heart ☐ Line Noise ☐ Chan Noise ☐ Other ☐ ?

4.




☐ Brain ☐ Muscle ☐ Eye ☐ Heart ☐ Line Noise ☐ Chan Noise ☐ Other ☐ ?

5.



☐ Brain ☐ Muscle ☐ Eye ☐ Heart ☐ Line Noise ☐ Chan Noise ☐ Other ☐ ?

6.



☐ Brain ☐ Muscle ☐ Eye ☐ Heart ☐ Line Noise ☐ Chan Noise ☐ Other ☐ ?

Exercise 3: P300 Evoked Potentials


The file `BI5_segments_HTS.mat` contains data from an auditory P300 BCI experiment. Study participants were asked to concentrate on a target tone, which was higher than the non-target tone. The multidimensional matrix `segments` stores the following information:

- First dimension: Channels (Fz, Cz and Pz, you can look that up in “ch_selection”).
- Second dimension: Time course (77 time points, Sample rate: 64 Hz).
- Third dimension: Trials (2200).

Use the `classlabels` vector to identify non-target trials (class 1) and target trials (class 2). Store them in separate matrices and calculate their means.

Tasks:

- Compute the EPs by averaging the trials for each channel and modality (non-target and target) separately.
- Plot a figure for each channel with the two modalities (`hold on`).
- Describe the differences between the two modalities (P300).
- Are the differences the same for all channels?



EASYCAP

EEG Recording Caps and Related Products


EASYCAP GmbH
Steingrabenstrasse 14
DE-82211 Herrsching
Germany
Delivery Address: Walchstadter Strasse 8a, 82266 Inning a. A.

Tel +49 (0) 8143 992-977
Fax +49 (0) 8143 992-979
www.easycap.de
info@easycap.de

Montage No. 15

See coordinates measured¹ from a real head in "M15-ext-ThetaPhi.xls".

Electrode Names:



Electrode positions:

| | | |
|---------|----------|------------|
| 1 Fp1 | 51 F7 | 101 C6 |
| 2 Fpz | 52 FT9 | 102 CCP8 |
| 3 Fp2 | 53 FAF5 | 103 TP8 |
| 4 AFp1 | 54 F5 | 104 P10 |
| 5 AFp2 | 55 FFC7 | 105 CFC4 |
| 6 AF3 | 56 FT7 | 106 C4 |
| 7 AFz | 57 F3 | 107 CCP6 |
| 8 AF4 | 58 FFC5 | 108 CP6 |
| 9 FAF1 | 59 FC5 | 109 PCP8 |
| 10 FAF2 | 60 CFC7 | 110 P8 |
| 11 F1 | 61 T7 | 111 PO10 |
| 12 Fz | 62 FFC3 | 112 CCP4 |
| 13 F2 | 63 FC3 | 113 CP4 |
| 14 FFC1 | 64 CFC5 | 114 PCP6 |
| 15 FFC2 | 65 C5 | 115 P6 |
| 16 FC1 | 66 CCP7 | 116 PO8 |
| 17 FCz | 67 TP7 | 117 OPO10 |
| 18 FC2 | 68 P9 | 118 PCP4 |
| 19 CFC1 | 69 CFC3 | 119 P4 |
| 20 CFC2 | 70 C3 | 120 PPO6 |
| 21 C1 | 71 CCP5 | |
| 22 Cz | 72 CP5 | 121 EMG1SL |
| 23 C2 | 73 PCP7 | 122 EMG1IL |
| 24 CCP1 | 74 P7 | 123 EMG2SL |
| 25 CCP2 | 75 P09 | 124 EMG2IL |
| 26 CP1 | 76 CCP3 | 125 EMG1SR |
| 27 CPz | 77 CP3 | |
| 28 CP2 | 78 PCP5 | |
| 29 PCP1 | 79 P5 | |
| 30 PCP2 | 80 PO7 | |
| 31 P1 | 81 OPO9 | |
| 32 Pz | 82 PCP3 | |
| 33 P2 | 83 P3 | |
| 34 PPO1 | 84 PPO5 | |
| 35 PPO2 | 85 AF8 | |
| 36 PO3 | 86 F10 | |
| 37 POz | 87 F8 | |
| 38 PO4 | 88 FT10 | |
| 39 OPO1 | 89 FAF6 | |
| 40 OPO2 | 90 F6 | |
| 41 O1 | 91 FFC8 | |
| 42 Oz | 92 FT8 | |
| 43 O2 | 93 F4 | |
| 44 OI1 | 94 FFC6 | |
| 45 OI2 | 95 FC6 | |
| 46 I1 | 96 CFC8 | |
| 47 Iz | 97 T8 | |
| 48 I2 | 98 FFC4 | |
| 49 AF7 | 99 FC4 | |
| 50 F9 | 100 CFC6 | |